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Purity of produced water: hig-
her than 200.000 ohm/cm for
water whose salt content does
not exceed 350 ppm.

IDROTECNICA produces devices
of higher capacity than those
listed.

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Automatic dual-column water demineralizers

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The Model 2CA water demineralizers consist of two separate packed-bed columns one containing a cationic, the other an anionic exchange resin with high ion exchange capacity. Unlike other very complicated units, the Model 2CA demineralizers are built to be simple, functional and reliable with tested, long-lasting components needing minimal maintenance. All components are made of acid-resistant materials. The ion exchange resins are packed in vinyl ester columns resistant to water pressure. An electronic conductivity meter, equipped with a variable threshold alarm, ensures constant control of the quality of the treated water, and indicates when resin regeneration needs to be performed. On request, the demineralizers may be equipped with a device to automatically stop delivery, or start discharge, of the water of insufficient purity.

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- Feed
- Cosmetics
- Chemistry
- Glass industry
- Galvanization
- Laboratory
- Hospital

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Resin regeneration is initiated by simply pressing a button, and occurs in fully automated fashion. On request, regeneration starts automatically when water purity falls below a preset value.

All demineralizers are equipped with two plastic containers for the regenerants. Regenerants are automatically withdrawn from the containers once they are placed on the floor. In case of power failure, resin regeneration may be performed manually. Idrotecnica shall analyze the water you have to treat free of charge to determine the demineralizer capacity that best suits your needs.