Demineralizzatori acqua a resine intercambiabile Water demineralizers with interchangeable resins

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IDROTECNICA

SERIE M

Semplice ed efficace.
Simple and efficient.

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The water demineralizers of the M Serie represent the simplest and most practical way to obtain cheaply demineralized water at high quality. The units of this serie have a moderate encumbrance, are very economical and can be used easily.

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The M Serie demineralizers are composed of a stiffened non-toxic polypropylene column, resistant to a pressure of 10 bar (M4: 1,5 bar).

The column contains selected cationic and anionic resins, having an high exchange power and treated and mixed in order to obtain water with a high purity.

According to the versions there are two different conductivity meters:

- a two led (red – green) conductivity meter to show when the resins are saturated, self-fed through a battery;
- a digital conductivity meter (scale 0–19,99 µS) with leds showing resins saturation for a continual checking of water purity on outlet, fed through a socket transformer 230/50Hz. All versions offer a kit of pressure-resistant connections and pipes with a sampling valve.

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You can effect the resins re-charge in two different ways:

- by replacing the resin charge, contained inside the column, with a new one supplied in a little sack (a very simple and quick operation);
- by replacing the column with a new one;
- the used resins and columns can be sent to our laboratory for a careful regeneration;
- M4 has a disposable cartridge.

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For the M Serie are available No.6 units (M4, M 9, M 12, M 20, M 30, M 50) in No.3 versions:

- endowed with connections kit, without conductivity meter
- endowed with two-led conductivity meter
- endowed with digital conductivity meter (not for M4).

In order to simplify the identification of the various models and spare parts, we created an »intelligent« code system by four numbers preceded by letter »M«. The first two numbers identify the column type (04, 09, 12, 20, 30, 50) the last two show the version and the spare part:

M-02 demineralizer without conductivity meter, only connections kit and sampling valve
M-04 demineralizer with two-led conductivity meter, connections kit and sampling valve
M-06 demineralizer with digital conductivity meter, connections kit and sampling valve
M-01 spare column (container + resin)
M-07 regenerated resins charge (in sack)
M-08 resin column regeneration service

Some examples:
Code **M0904**
= demineralizer M 9 with two-led conductivity meter
Code **M3006**
= demineralizer M 30 with digital conductivity meter,
Code **M5001**
= spare column for M 50